

620100 Surgical Hood Tuck

Surgical Hood

Product details

Size :	One size
Descriptive feature :	Standard, Tie Band
Colour :	Green
Sterile :	Non-sterile
Brand :	BARRIER

Images



Delivered items

620100-11

Regulatory released: Algeria, Austria, Belgium, Bulgaria, Côte d'Ivoire, Cyprus, Czechia, Denmark, Estonia, Finland, France, French Polynesia, Germany, Greece, Guadeloupe, Hong Kong, Hungary, Italy, Latvia, Lithuania, Luxembourg, Malaysia, Martinique, Mauritius, Moldova (the Republic of), Monaco, Morocco, Netherlands, New Caledonia, Norway, Oman, Poland, Portugal, Qatar, Romania, Russian Federation, Saudi Arabia, Singapore, Slovakia, Spain, Sweden, Switzerland, Turkey, United Arab Emirates, United Kingdom of Great Britain and Northern Ireland

Country of origin: China

Shelf life: 5 years

Sterilization method: Non-sterile

Packing information: First packaging layer is a cardboard dispenser box. Second layer is a corrugated board transport box.

Is suitable for Tray: Yes

Packing level	Quantity	GS1 code	WxLxH (mm)	Vol (dm ³)	Weight gross/net (kg)
Dispenser box	100	7332430183910	215x263x150		
Transport box	500	7310791008028	280x760x230	48,9	4,9 / 3,5

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2021 Mölnlycke Health Care AB. All rights reserved.



Packing level	Quantity	GS1 code	WxLxH (mm)	Vol (dm3)	Weight gross/net (kg)
Pallet	14000	7332430548726	800x1200x1754		

620100-40

Regulatory released: Algeria, Austria, Belgium, Bulgaria, Cyprus, Czechia, Denmark, Estonia, Finland, France, French Polynesia, Germany, Greece, Guadeloupe, Hungary, Italy, Latvia, Lithuania, Luxembourg, Martinique, Mauritius, Moldova (the Republic of), Monaco, Netherlands, New Caledonia, Norway, Oman, Poland, Portugal, Qatar, Romania, Saudi Arabia, Singapore, Slovakia, Slovenia, Spain, Sweden, Switzerland, Turkey, United Arab Emirates, United Kingdom of Great Britain and Northern Ireland

Country of origin: Morocco

Shelf life: 5 years

Sterilization method: Non-sterile

Packing information: First packaging layer is a cardboard dispenser box. Second layer is a corrugated board transport box.

Is suitable for Tray: No

Packing level	Quantity	GS1 code	WxLxH (mm)	Vol (dm3)	Weight gross/net (kg)
Dispenser box	120	7323190221219	145x252x210		
Transport box	960	7323190221202	265x593x439	69,0	6,3 / 4,8
Pallet	23040	7323190221196	800x1200x1906		176,7 / 115,2

620100-41

Regulatory released: Austria, Belgium, Bulgaria, Cyprus, Czechia, Denmark, Estonia, Finland, France, French Polynesia, Germany, Greece, Hong Kong, Hungary, Italy, Latvia, Lithuania, Luxembourg, Martinique, Moldova (the Republic of), Monaco, Netherlands, Norway, Oman, Poland, Portugal, Qatar, Romania, Saudi Arabia, Singapore, Slovakia, Spain, Sweden, Switzerland, United Arab Emirates, United Kingdom of Great Britain and Northern Ireland

Country of origin: China

Shelf life: 5 years

Sterilization method: Non-sterile

Packing information: First packaging layer is a cardboard dispenser box. Second layer is a corrugated board transport box.

Is suitable for Tray: No

Packing level	Quantity	GS1 code	WxLxH (mm)	Vol (dm3)	Weight gross/net (kg)
Dispenser box	120	7323190275199	145x252x210		0,8 / 0,7
Transport box	960	7323190275182	265x613x439	71,3	7,5 / 5,4
Pallet	23040	7323190275175	800x1200x1906		205,0 / 129,1

Material

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2021 Mölnlycke Health Care AB. All rights reserved.



Animal tissues :	No
Human blood derivatives :	No
Natural rubber latex :	No
Medicinal substances :	No
Phthalates :	No
Polyvinyl chloride :	No

Product Composition Headwear

Product Component	Composition
Main material	Viscose nonwoven
Sewing thread	Polyester

Technical

Dimension

Dimension text	Dimension value
Length straps	62 cm
Length side piece	28 cm
Face length	17 cm
Face width	9,5 cm
Width tie-band	4 cm
Length top piece	25 cm
Width top piece	15 cm
Height forehead	4 cm

Disposal instructions

Non-hazardous waste used BARRIER products and sterility barriers should, in the majority of cases, be classified as non-hazardous waste. They contain high amounts of energy and are well suited for incineration. BARRIER products do not contain any hazardous substances that can leach out if the products are land filled. Transport boxes are designed to fit existing recovery systems. The new BARRIER packaging system complies with the Packaging Waste Directive of the European Union.

Storage instructions

Mölnlycke Health Care recommends that BARRIER products are stored under normal storage conditions. All layers of packaging should be kept intact until access to the underlying layers is needed. Storage facilities for products only protected by the sterility barrier should be kept under conditions where low level of particulate air contamination prevail, so that it would not constitute a risk to the patient when the package is opened and the product is used.

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2021 Mölnlycke Health Care AB. All rights reserved.



Classifications

Regulation type(s)	MDD Class I	Locally Regulated	Unregulated
Intended use MDD :	The intended use of headwear is infection control by reducing the risk of transferring bacteria, hair or particles into the surgical wound.		
MDD Classification Rule :	1		
Conformity Annexes :	VII		
Measuring Function :	No		
Regulatory released :	Austria, Belgium, Bulgaria, Cyprus, Czechia, Denmark, Estonia, Finland, France, Germany, Greece, Hungary, Italy, Latvia, Lithuania, Luxembourg, Martinique, Monaco, Netherlands, New Caledonia, Norway, Poland, Portugal, Romania, Slovakia, Slovenia, Spain, Sweden, Switzerland, United Kingdom of Great Britain and Northern Ireland	Algeria, Guadeloupe, Malaysia, Moldova (the Republic of), Morocco, Oman, Qatar, Russian Federation, Saudi Arabia, Turkey, United Arab Emirates	Côte d'Ivoire, French Polynesia, Hong Kong, Mauritius, Singapore

Applied standards : The standards presented below is a selection of the most essential standards that are adhered to.

EN 1041, EN 62366, EN ISO 9001, EN ISO 13485, EN ISO 10993-1, EN ISO 10993-5, EN ISO 15223-1, ISO 15223-2, EN ISO 10993-10, ISO 14001

Removable label

No

GMDN Code (Global Medical Device Nomenclature)

32297 Surgical cap, single-use

Reusability

Single use

UNSPSC

42131611 Surgeon caps or hoods

Commodity Code

6505009090 Caps and hoods

Find out more at www.molnlycke.com

